

Research Product Order Form

-Complete the form below for <u>each separate</u> delivery date and email to <u>Info@BioSharingNetwork.com</u>. -Please submit orders at least <u>5 business days</u> prior to the time products will be shipped or picked up.

-Ship/pick-up times are **4:00pm** unless otherwise stated in the form below.

-Orders submitted 3 or 4 days before ship/pickup time will incur a **\$25 stat fee**.

-Orders submitted 2 days before ship/pickup time will incur a **\$40 stat fee**.

-Orders received 1 day before ship/pickup date **<u>cannot</u>** be filled.

-For buffy coat and platelet orders on days with low predicted whole blood draws, you may be asked to choose a different ship date or different product.

| Re | equesting Company: | | | | | | | | |
|---|---|---|------------------|----------------------------|----------|--|--|--|--|
| Primary C | Contact and/or Lab Name: | | | | | | | | |
| Pe | erson Placing Order: | | | | | | | | |
| | Delivery Date: | | | | | | | | |
| Product Use: Research Further Manufacturing Other | | | | | | | | | |
| Catalog Number | Description | | Product Price | Donor Demo- graphics | Quantity | | | | |
| BUFFY COATS - Whole blood-derived, manufactured <8 hours of collection; 1-2.5x10 ⁹ leukocytes (30-80% PMNs) in ~50 mL plasma | | | | | | | | | |
| BC-8L | Buffy Coat; available at 4:00pm day of collection | | | \$100 | \$10 | | | | |
| BC-8E | Buffy Coat; available at 10:00am day after collection | | | \$100 | \$10 | | | | |
| LEUKOREDUCTION SYSTEM (LRS) CHAMBERS - Byproduct of platelet apheresis; concentrated PBMCs (1-2.5x10 ⁹ PBMCs, <10% PMNs) in ~10 mL blood; "buffy cone" | | | | | | | | | |
| LRSC-01 | LRS chamber, sterilely sealed; available at 4:00pm day of collection | | \$80 | \$10 | | | | | |
| LRSC-02 | LRS chamber, sterilely sealed; available at 12:00pm day after collection | | | \$65 | Included | | | | |
| LRSC-03 | LRS chamber; drained, cells suspended in buffer; >1.2x10 ⁹ PBMC ** | | | \$90 | Included | | | | |
| ISOLATED PBMCs - Isolated by Ficoll from LRS chambers the day after collection ** Available only Tuesday – Friday after 12:00pm. Cell counts cannot be guaranteed when specific donor demographics are requested. | | | | | | | | | |
| PBMC-01 | PBMCs suspended in buffer, >0.8x10 ⁹ PBMCs ** | | \$200 | Included | | | | | |
| PBMC-02 | PBMCs, cryopreserved; 10 ⁸ /v | al | | \$125 | Included | | | | |
| PLATELETS - Whole blood-derived platelet rich plasma (PRP); CPD anticoagulant | | | | | | | | | |
| PRP-8S | Single-donor "random" PRP; | avg of ~6x10 ¹⁰ platelets in 40-60 mL | plasma | \$100 | \$10 | | | | |
| PRP-8C | Platelet count provided on si | ngle-donor PRP unit | | \$25 | N/A | | | | |
| PRP-24P | Pooled PRP from 5-6 donors; | avg of ~3x10 ¹¹ platelets in 250-350 n | nL plasma | \$200 | N/A | | | | |
| PLASMA - Whole blood-derived, 250-300 mL/unit; CPD anticoagulant | | | | | | | | | |
| RPFF | Plasma isolated from whole b | lood and frozen <24 hours after co | llection | \$60 | \$10 | | | | |
| RED BLOOD CELL (RBC) AND WHOLE BLOOD UNITS - (subject to availability) | | | | | | | | | |
| RBC-01 | Packed RBC unit, fresh, Adso | anticoagulant; available day after | collection | \$280 | \$10 | | | | |
| RBC-02 | Packed RBC unit, expired; Ac | lsol anticoagulant | | \$25 | \$10 | | | | |
| WB-01 | Whole blood unit, fresh; CPD | anticoagulant; available day of cc | llection | \$550 | \$10 | | | | |

| DTA-01 | Single 6 mL purple-top tube containing whole blood in EDTA anticoagulant | \$20 | \$10 |
|----------|---|-------------|-------------|
| EDTA-02 | Two 6 mL EDTA purple-top tubes from same donor | | \$10 |
| SRT-01 | Single 6 mL red-top tube containing clotted whole blood | | \$10 |
| SRT-02 | Two 6 mL red top tubes from the <u>same donor</u> | | \$10 |
| SRT-03 | Pooled serum only from 2 red top tubes from the <u>same donor</u> (4-5 mL) | | Included |
| | JBE SERUM, MULTI-DONOR - onors/L of serum, pooled 72-96 hours from collection, verified negative infectiou | s disease t | est results |
| PTS-01 | 3.5 L Pooled tube serum, not sterile, frozen | | N/A |
| PTS-02 | 1.0 L Pooled tube serum, not sterile, frozen | | N/A |
| PTS-03 | 1.0 L Pooled tube serum, not sterile, refrigerated (made on-demand) | | N/A |
| PTS-04 | 1.0 L Pooled tube serum, filter sterilized, refrigerated (made on-demand) | | N/A |
| PTS-05 | 1.0 L Pooled tube serum, filter sterilized, frozen | | N/A |
| PTS-06 | 500 mL Pooled tube serum, filter sterilized, frozen | | N/A |
| PTS-07 | 250 mL Pooled tube serum, filter sterilized, frozen | | N/A |
| PTS-08 | 50 mL Pooled tube serum, filter sterilized, frozen | | N/A |
| DONOR SE | LECTION AND RECRUITMENT - | | |
| Sel-01 | For the selection of an already-collected product based on specific donor attributes Age: Gender: Race: CMV: Other: Other: CMV: CMV: | Quote | Included |
| Rec-01 | For the recruitment of a donor for a product collection based on specific demographics Age: Gender: Race: CMV: Other: Other: CMV: | Quote | Included |

Shipping/Pick-up Information Product(s) will be picked up at BSN Center. Phone number to call/text when order is ready for pick up: For pick-up orders, please fill out all fields below except FedEx Account Number and Delivery Time. Note: A shipping charge will be added if no FedEx account number is provided FedEx Account Number: Shipping Temperature: Ambient Refrigerated Frozen **Delivery Time:** 10:30am 8:30am Ship-to Company: Shipping Address: City: State: Zip: Contact Person: Phone Number: Billing Information - PO or Credit Card (Visa/MasterCard only) PO Number: CC Name: CC Number: Exp. Date: CVV Code: